



DEFENCE



DÉFENSE

# **Electrospray Mass Spectrometry of Chemical Warfare Agents, Degradation Products and Related Compounds**

Paul A. D'Agostino, Claude L. Chenier and James R. Hancock  
Defence Research Establishment Suffield

**Defence R&D Canada**

Technical Report

DRES TR 2002-028

March 2002

**DISTRIBUTION STATEMENT A**  
Approved for Public Release  
Distribution Unlimited



National  
Defence

Défense  
nationale

**Canada**

**20020520 153**

# **Electrospray Mass Spectrometry of Chemical Warfare Agents, Degradation Products and Related Compounds**

Paul A. D'Agostino  
Defence Research Establishment Suffield

Claude L. Chenier  
Defence Research Establishment Suffield

James R. Hancock  
Defence Research Establishment Suffield

**Defence Research Establishment Suffield**

Technical Report  
DRES TR 2002-028  
March 2002

Authors

---

Paul A. D'Agostino, Claude L. Chenier and James R. Hancock

Approved by

---

Cam Boulet

H/CBDS

Approved for release by

---

Robin Clewley

DRP Chairperson

## Abstract

---

Over the past several years Defence Research Establishment Suffield (DRES) developed LC-ESI-MS methods for the analysis of a variety of samples of chemical defence interest. These samples included hydrolysed munitions grade mustard, degraded O-ethyl S-[2-(diisopropylamino)ethyl] methylphosphonothiolate (VX), snow contaminated with sarin (GB), soil contaminated with sarin (GB) and soman (GD), munitions grade tabun (GA), failed synthetic products and synthetic standards. During these investigations DRES acquired ESI-MS data for sixty chemical warfare agents, degradation products and related compounds. Data were obtained with a resolution of 5000 (50% valley definition) in the continuum mode at several sampling cone voltages. At lower sampling cone voltages the mass spectra were dominated by protonated, sodiated and protonated acetonitrile adducts and/or their dimers that could be used to confirm the molecular mass of each compound. Structural information was obtained by inducing product ion formation in the ESI interface at higher sampling cone voltages. Representative data obtained at both lower and higher sampling cone voltages were selected for entry into the DRES ESI-MS Database as part of a collaboration between Canada (DRES) and The Netherlands (TNO Prins Maurits Laboratory).

## Résumé

---

Durant ces dernières années, le Centre de recherches pour la défense, Suffield (CRDS) a mis au point des méthodes CL-ENI-SM (chromatographie liquide – ionisation par électronébulisation – spectrographie de masse) pour analyser une variété d'échantillons d'intérêt pour la défense chimique. Ces échantillons comprenaient de l'ypérite à munitions hydrolysée, du méthylphosphonothiolate de O-éthyle et de S (diisopropylamino-2 éthyle) (VX) dégradé, de la neige contaminée par du sarin (GB), du sol contaminé par du sarin (GB) et du soman (Gd), du tabun à munitions (GA), des produits synthétiques incomplets et des normes synthétiques. Durant ces études, le CRDS a acquis des données SM-ENI pour soixante agents de guerre chimique, produits de dégradation et composés reliés. Les données ont été obtenues avec une résolution de 5000 (à un niveau équivalent à 50% du pic chromatographique) en mode continu des tensions des entonnoirs de prélèvement. À des tensions plus faibles, les spectres de masse étaient dominés par des adduits protonés, de sodium et acétonitriles protonés et/ou des dimères qui pourraient être utilisés pour confirmer la masse moléculaire de chaque composé. L'information structurale a été obtenue en induisant la formation des ions du produit dans l'interface ENI à des tensions plus hautes des entonnoirs de prélèvement. Des données représentatives obtenues avec des tensions d'entonnoirs de prélèvement à la fois plus hautes et plus faibles ont été sélectionnées pour être entrées dans la base de données ENI-SM du CRDS durant les travaux de collaboration entre le Canada (CRDS) et les Pays-Bas (laboratoire Prins Maurits TNO).

This page intentionally left blank.

## Executive summary

---

**Introduction:** The Canadian Forces (CF) and its allies may be called on to perform peacekeeping or battlefield operations in regions of the world where there is a significant threat of chemical/biological (CB) warfare agent use. To operate effectively in these theatres they must be able to identify the CB agent used. Mass spectrometry (MS) is a powerful analytical technique for the identification of both known and unknown compounds and DRES is currently investigating this instrumental technique in fulfillment of CF CB detection and identification requirements.

**Results:** Over the past several years DRES developed liquid chromatography-electrospray-mass spectrometry (LC-ESI-MS) methods for the analysis of a variety of contaminated samples. During the course of these investigations DRES acquired ESI-MS data for sixty chemical warfare agents, degradation products and related compounds. Representative data were selected for entry into the DRES ESI-MS Database as part of a collaboration between Canada (DRES) and The Netherlands (TNO Prins Maurits Laboratory).

**Significance:** The CF and its allies may be deployed in regions of the world where there is a significant threat of chemical/biological warfare agent use. Gas chromatography-mass spectrometry, while suitable for the direct analysis of organophosphorus chemical warfare agent in organic extracts, is usually not preferred for the direct analysis of aqueous samples or extracts since these samples normally require additional sample handling steps and derivatization prior to analysis. LC-ESI-MS is being used increasingly, as electrospray mass spectrometric data may be used to directly identify chemical warfare agents, degradation products and related compounds in collected aqueous samples or extracts. Identification of the CB agent is of importance since the results of such analyses would contribute to the development of strategic and political positions regarding future military operations and would facilitate the dissemination of technical advice to in-theatre field commanders and medical personnel.

**Future Plans:** The mass spectra reported in this database would be valuable for the identification of chemical warfare agents, degradation products and related compounds samples collected by the military or in support of Chemical Weapons Convention inspections. The database will be updated as part of the on-going collaboration between Canada and The Netherlands.

D'Agostino, P.A., Chenier, C.L., Hancock, J.R. 2002. Electrospray Mass Spectrometry of Chemical Warfare Agents, Degradation Products and Related Compounds. DRES TR 2002-028. Defence Research Establishment Suffield.

## Sommaire

---

**Introduction:** Les Forces canadiennes (FC) et leurs alliés peuvent être appelés à exécuter des opérations de maintien de la paix et des opérations sur les champs de bataille dans les régions du monde où il existe une menace sérieuse d'utilisation d'agents chimiques et biologiques (CB). Pour opérer de manière efficace dans ces théâtres, les FC doivent être en mesure d'identifier les agents CB utilisés. La spectrométrie de masse (SM) est une technique analytique d'identification des composés à la fois connus et inconnus. Le Centre de recherches pour la défense, Suffield (CRDS) est en train d'étudier cette technique instrumentale dans le but de répondre aux exigences des FC de détection et d'identification CB.

**Résultats :** Ces dernières années, le CRDS a mis au point des méthodes CL-ENI-SM (chromatographie liquide – ionisation par électronébulisation – spectrographie de masse) pour analyser une variété d'échantillons contaminés. Durant le cours des recherches, le CRDS a acquis des données SM-ENI pour soixante agents chimiques de guerre, produits de dégradation et composés reliés. Les données représentatives ont été sélectionnées pour être entrées dans la base de données SM -ENI du CRDS durant les travaux de collaboration entre le Canada (CRDS) et les Pays-Bas (laboratoire Prins Maurits TNO).

**La portée des résultats :** Les FC et leurs alliés peuvent être appelés à se déployer dans les régions du monde où il existe une menace sérieuse d'utilisation des agents de guerre chimique et biologique. La spectrométrie de masse de chromatographie des gaz convient pour l'analyse directe des agents de guerre chimique organophosphorés mais elle n'est cependant pas la meilleure méthode pour l'analyse directe des échantillons aqueux ou des extraits, étant donné que ces échantillons requièrent normalement des étapes de manipulation d'échantillons additionnelles et une dérivation antérieure à l'analyse. La méthode CL-ENI-SM est de plus en plus utilisée, les données SM-ENI pouvant être utilisées pour identifier les agents de guerre chimique, les produits de dégradation et les composés reliés directement à partir d'échantillons aqueux ou d'extraits. L'identification des agents CB est de grande importance puisque les résultats de telles analyses contribueront au développement des positions stratégiques et politiques des futures opérations militaires et faciliteront la dissémination des conseils techniques aux commandants sur le champ des théâtres d'opérations et au personnel médical.

**Les buts futurs :** Les spectres de masse introduits dans les bases de données seront précieux pour l'identification des agents de guerre chimique, les produits de dégradation et les échantillons composés reliés cueillis par le militaire ou en soutien aux inspections de la Convention sur les armes chimiques. Étant donné la collaboration continue entre le Canada et les Pays-Bas, la base de donnée continuera d'être mise à jour.

D'Agostino, P.A., Chenier, C.L., Hancock, J.R. 2002. Electrospray Mass Spectrometry of Chemical Warfare Agents, Degradation Products and Related Compounds. DRES TR 2002-028. Defence Research Establishment Suffield.

## Table of contents

---

Abstract .....	i
Résumé .....	i
Executive summary .....	iii
Sommaire .....	iv
Table of contents .....	v
List of figures .....	vi
List of tables .....	viii
Introduction .....	1
Experimental .....	3
Results and Discussion .....	4
LC-ESI-MS separations .....	4
LC-ESI-MS sensitivity .....	7
ESI-MS mass spectral differences .....	7
Conclusions .....	72
References .....	73



## List of figures

---

Figure 1.	Diisopropylamine (101) .....	12
Figure 2.	2,2'- Thiodiethanol (Thiodiglycol) (122).....	13
Figure 3.	Ethyl methylphosphonic acid (124) .....	14
Figure 4.	Isopropyl methylphosphonic acid (138).....	15
Figure 5.	Ethyl ethylphosphonic acid (138) .....	16
Figure 6.	Isopropyl methylphosphonofluoridate (Sarin, GB) (140) .....	17
Figure 7.	Diethyl methylphosphonate (152).....	18
Figure 8.	Isopropyl methyl methylphosphonate (152) .....	19
Figure 9.	Ethyl dimethylphosphoramidic acid (153).....	20
Figure 10.	Tetramethyl phosphorodiamidic cyanide (161) .....	21
Figure 11.	O-Ethyl N,N-dimethylphosphoroamidocyanidate (Tabun, GA) (162) .....	22
Figure 12.	Cyclohexyl methylphosphonic acid (178).....	23
Figure 13.	Hexamethylphosphorotriamide (179) .....	24
Figure 14.	Cyclohexyl methylphosphonofluoridate (GF) (180).....	25
Figure 15.	Pinacolyl methylphosphonic acid (180).....	26
Figure 16.	Diisopropyl methylphosphonate (180).....	27
Figure 17.	Ethyl tetramethylphosphorodiamidate (180).....	28
Figure 18.	Diethyl dimethylphosphoramidate (181) .....	29
Figure 19.	3,6-Dithia-1,8-octanediol (182).....	30
Figure 20.	Triethyl phosphate (182) .....	31
Figure 21.	Pinacolyl methylphosphonofluoridate (Soman, GD) (182) .....	32
Figure 22.	8-Chloro-6-oxa-3-thia-1-octanol (184) .....	33
Figure 23.	Diisopropyl phosphorofluoridate (DFP) (184).....	34

Figure 24.	Tributylamine (185) .....	35
Figure 25.	2-(Diisopropylamino)ethyl ethyl sulfide (189) .....	36
Figure 26.	Diisopropyl ethylphosphonate (194) .....	37
Figure 27.	Ethyl isopropyl dimethylphosphoramidate (195) .....	38
Figure 28.	Diethyl isopropyl phosphate (196) .....	39
Figure 29.	Pentyl dibutylamine (199) .....	40
Figure 30.	2,2',2''-Trichlorotriethylamine (HN3) (203) .....	41
Figure 31.	1,7-Dioxa-4,10-dithiacyclodecane (208) .....	42
Figure 32.	Diisopropyl dimethylphosphoramidate (209) .....	43
Figure 33.	Diisopropyl ethyl phosphate (210) .....	44
Figure 34.	Diisopropyl chloromethylphosphonate (214) .....	45
Figure 35.	Triisopropyl phosphate (224) .....	46
Figure 36.	N,N'-Dicyclohexylurea (224) .....	47
Figure 37.	6-Oxa-3,9-dithia-1,11-undecanediol (226) .....	48
Figure 38.	Diethyl dimethylpyrophosphonate (230) .....	49
Figure 39.	S-[2-(Diisopropylamino) ethyl]methylphosphonothiolate (239) .....	50
Figure 40.	N,N'-Dicyclohexylthiourea (240) .....	51
Figure 41.	3,6,9-Trithia-1,11-undecanediol (242) .....	52
Figure 42.	Diisopropyl dichloromethylphosphonate (248) .....	53
Figure 43.	Ethyl phosphoric tetramethylphosphorodiamidic anhydride (260) .....	54
Figure 44.	Diethyl phosphoric ethyl hydrogen phosphoric anhydride (262) .....	55
Figure 45.	O-Ethyl S-[2-(diisopropylamino)ethyl] methylphosphonothiolate (VX) (267) .....	56
Figure 46.	6-Oxa-3,9,12-trithia-1,4-tetradecanediol (286) .....	57
Figure 47.	Bis(tetramethylphosphorodiamidic) anhydride (286) .....	58
Figure 48.	Ethyl dimethylphosphoramidic tetramethylphosphoamidic anhydride (287) .....	59

Figure 49.	Bis[2-(diisopropylamino)ethyl] sulfide (288) .....	60
Figure 50.	Bis(ethyl dimethylphosphoramidic) anhydride (288) .....	61
Figure 51.	Diethyl phosphoric tetramethylphosphorodiamidic anhydride(288) .....	62
Figure 52.	Diethyl phosphoric ethyl dimethylphosphoramidic anhydride (289).....	63
Figure 53.	Bis(diethyl phosphoric) anhydride (TEPP) (290) .....	64
Figure 54.	Diethyl phosphoric ethyl isopropyl phosphoric anhydride (304).....	65
Figure 55.	Bis[2-diisopropylamino)ethyl] disulfide (320) .....	66
Figure 56.	O-Ethyl-S-[5-(diisopropylamino)-3-thiapentyl] methylphosphonothiolate (327).....	67
Figure 57.	Diisopropyl phosphoric ethyl isopropyl phosphoric anhydride (332).....	68
Figure 58.	Bis(diisopropylphosphoric) anhydride (346) .....	69
Figure 59.	Octamethyltetramidotriphosphoric acid (or isomer) (366) .....	70
Figure 60.	O-Ethyl-S[8-(diisopropylamino)-3,6-dithiaoctyl] methylphosphonothiolate (387)...	71

## List of tables

---

Table 1.	Compounds identified in contaminated snow by LC-ESI-MS (23).....	5
Table 2.	Principal compounds identified in munitions-grade tabun sample (24). .....	6
Table 3.	Lower mass product ions observed during chemical warfare agent analysis by LC-ESI-MS. ....	9
Table 4.	Typical high resolution ESI-MS data acquired for compounds identified in the contaminated snow sample. (Sampling cone voltage: 20 volts) (23).....	10
Table 5.	High resolution ESI-MS data for two compounds identified in a failed tabun synthesis sample (Sampling cone voltage: 40 volts) (26).....	11

## Introduction

---

Canada, one of more than 140 countries to ratify the Chemical Weapons Convention, complies with this treaty by maintaining only small stocks of chemical warfare agents at Defence Research Establishment Suffield (DRES) for defensive research and development. DRES maintains the Canadian National Single Small-Scale Facility and provides the Canadian Forces and other National interests with a laboratory based identification capability for chemical warfare agents, degradation products and related compounds. Analyses of this type requires the use of sensitive, specific analytical methods, particularly when unambiguous proof is required for the presence of chemical warfare agents. These analytical demands are being actively addressed by DRES through the development and application of new analytical methods for the detection and identification of chemical warfare agents in a variety of samples.

Gas chromatography (GC) has been used extensively for the separation and identification of chemical warfare agents, with gas chromatography-mass spectrometry (GC-MS) being used frequently for the characterization of these compounds (1,2). GC-MS, while suitable for the direct analysis of organophosphorus chemical warfare agent in organic extracts, is usually not preferred for the direct analysis of aqueous samples or extracts since these samples normally require additional sample handling steps and derivatization prior to analysis. LC-ESI-MS is being used increasingly, as electrospray mass spectrometric data may be used to directly identify chemical warfare agents, degradation products and related compounds in collected aqueous samples or extracts. Increasingly researchers have developed atmospheric pressure ionization (e.g., electrospray (ESI), ionspray and atmospheric pressure chemical ionization) methods for the characterization of polar pesticides (3), organophosphate esters (4), and chemical warfare agents and/or their degradation products (5-24). These ionization modes have been interfaced to liquid chromatography (LC) and capillary electrophoresis (CE), with LC-MS (9-12, 14, 15, 18-24) and CE-MS (5, 13) methods being reported for the identification of lower volatility chemical warfare agent hydrolysis products. Recently, several LC-ESI-MS papers were published on the identification of organophosphorus chemical warfare agents and their hydrolysis products in aqueous (or snow) samples (14-17, 20, 23, 24) and aqueous extracts of contaminated soil samples (21, 22) during a single analysis. LC-ESI-MS was an attractive alternative to GC-MS for aqueous analyses, as both the organophosphorus chemical warfare agents and their hydrolysis products could be analysed directly without additional sample handling and derivatization.

DRES has analysed a variety of contaminated samples by LC-ESI-MS over the past several years using a time-of-flight mass spectrometer including:

- a) Hydrolysed munitions grade mustard
- b) Degraded O-ethyl S-[2-(diisopropylamino)ethyl] methylphosphonothiolate (VX)
- c) Snow contaminated with sarin (GB)
- d) Soil contaminated with sarin (GB) and soman (GD)

- e) Munitions grade tabun (GA)
- f) Failed synthesis of tabun (GA)
- g) Synthetic standards

The ESI-MS data acquired for the chemical warfare agents, degradation products and related compounds during these analyses were obtained with a resolution of 5000 (50% valley definition) in the continuum mode at several sampling cone voltages. At lower sampling cone voltages (typically 20 volts) the mass spectra were dominated by protonated, sodiated and protonated acetonitrile adducts and/or their dimers that could be used to confirm the molecular mass of each compound. Structural information was obtained by inducing product ion formation in the ESI interface at higher sampling cone voltages (typically 30 volts or more). Representative data obtained at both lower and higher sampling cone voltages for each compound were selected for entry into the DRES ESI-MS Database as part of a database collaboration between Canada (DRES) and The Netherlands (TNO Prins Maurits Laboratory) (25).

## Experimental

---

ESI-MS data were acquired using a Micromass LCT time-of-flight mass spectrometer equipped with the Z-spray electrospray interface. The electrospray capillary was operated at 3.2 kV with sampling cone voltages in the 20 to 70 volt range. Nitrogen desolvation gas was introduced into the interface (80 °C) at a flow rate of 480 L/h. Nitrogen nebulizer gas was introduced at a flow rate of 66 L/h. ESI-MS data were generally acquired from 700 to 70 Da (1 sec) in the continuum mode with a resolution of 5000 (50% valley definition).

LC separations were performed using an Applied Biosystems Model 140B dual syringe pump equipped with a Micro-Tech 150 mm × 0.32 mm I.D. C<sub>18</sub> (Zorbax, 5 μm) packed fused-silica capillary column and a Rheodyne 8125 injector with a 5 μL sample loop. The following solvent compositions were prepared for LC separation: Solvent A (0.1% trifluoroacetic acid in water) and Solvent B (0.1% trifluoroacetic acid in acetonitrile (ACN)/water, 95:5).

Chromatographic separations were performed using a 1% to 75%B gradient program over 30 (or 45) minutes or a 1% to 40%B gradient program over 30 minutes. In order to minimize dead volume effects and ensure reproducible mixing, the mobile phase was delivered at 200 μL/min and split prior to the injector such that the flow through the column was 16 μL/min.

## Results and Discussion

---

The DRES ESI-MS Database contains the mass spectra of 60 chemical warfare agents, degradation products and related compounds acquired during the analysis of chemical munitions fills, synthetic mixtures and contaminated environmental samples. The mass spectra, illustrated in Figures 1 to 60, have been compiled in order of increasing molecular mass and contain ESI-MS data acquired at two sampling cone voltages for each compound.

### LC-ESI-MS separations


Collected environmental and other samples typically contain multiple components and require chromatographic separation prior to MS characterization of each individual component. Capillary column GC-MS remains the method of choice for the analysis of many contaminated organic extracts, with electron impact ionization being used most frequently for the characterization of each compound (1,2). However this separation technique does not generally permit direct analysis of contaminated aqueous samples or extracts and requires a derivatization step prior to the analysis of many of the more polar chemical warfare agent degradation products. Increasingly, researchers have turned to LC-ESI-MS (14-17, 20, 23, 24) for the analysis of chemical warfare agents and their degradation products in aqueous samples, as this technique allows for the identification of these compounds in aqueous samples during a single analysis.

Several acetonitrile/water gradient programs were established to facilitate LC separation of the aqueous sample components in the samples studied. Flow rates, ranging from 10 to 20  $\mu\text{L}/\text{min}$  were evaluated, with 16  $\mu\text{L}/\text{min}$  being selected as a compromise, offering good chromatographic peak shape and component separation in a reasonable analysis time. Most samples contained multiple components, with the snow sample contaminated with sarin (Table 1) and the munitions grade tabun sample (Table 2) being typical. Upwards of a dozen sample components were readily identified and characterized in both examples by LC-ESI-MS (23, 24).

LC-ESI-MS offers several advantages over GC-MS for the analysis of chemical warfare agents, degradation products and related compounds. This analytical approach allows for the analysis of most chemical warfare agents (with the exception of mustard) and their more polar degradation (e.g., hydrolysis) products during a single analysis as opposed to the two analyses typically required by GC-MS. This advantage was recently demonstrated during the analysis of aqueous extracts of several soil samples contaminated with sarin and soman (21, 22). LC-ESI-MS was used to identify sarin (Figure 6) and soman (Figure 21) as well as their hydrolysis products, isopropyl methylphosphonic acid (Figure 4) and pinacolyl methylphosphonic acid (Figure 15). Related sample components were also readily identified in the presence of chemical warfare agents and degradation products, as illustrated by the number of compounds of similar structure characterized during LC-ESI-MS analysis of the contaminated snow sample and the munitions-grade tabun sample, respectively (Tables 1 and 2).

**Table 1.** Compounds identified in contaminated snow by LC-ESI-MS (23).

LC-ESI-MS Peak #	Compound Name	MW
1	Methylphosphonic acid	96
2	Methylphosphonofluoridic acid	98
3	Dimethylphosphinic acid	94
4	Isopropyl methylphosphonic acid	138
5	Isopropyl dimethylphosphinate	136
6	Sarin (Isopropyl methylphosphonofluoridate)	140
7	Isopropyl methyl methylphosphonate	152
8	Diisopropyl methylphosphonate	180
9	Pentyl dibutylamine	199
10	Tributylamine	185
11	Diisopropyl chloromethylphosphonate	214
12	Diisopropyl phosphorofluoridate	184
13	Triisopropyl phosphate	224
14	Diisopropyl dichloromethylphosphonate	248




LC-ESI-MS total-ion-current (600 to 85 Da) chromatogram obtained for the snow sample with a sampling cone voltage of 20 volts. Gradient program: 1% to 75%B in 30 minutes.



**Table 2.** Principal compounds identified in munitions-grade tabun sample (24).

LC-ESI-MS Peak #	Compound Name	MW
1	Ethyl dimethylphosphoramidic acid	153
2	Tetramethylphosphordiamidic cyanide	161
3	Ethyl dimethylphosphoramidocyanide (tabun)	162
4	Ethyl tetramethylphosphoramidate	180
5	Diethyl dimethylphosphoramidate	181
6	Triethyl phosphate	182
7	Diethyl phosphoric ethyl dimethylphosphoramidic anhydride	289
8	Ethyl isopropyl dimethylphosphoramidate	195
9	Diethyl isopropyl phosphate	196
10	Diethyl phosphoric isopropyl dimethylphosphoramidic anhydride (or isomer)	303
11	Diisopropyl ethyl dimethylphosphoramidate	209
12	Diethyl phosphoric ethyl isopropyl phosphoric anhydride	304
13	Diisopropyl ethyl phosphate	210
14	Diethyl phosphoric diisopropyl phosphoric anhydride (or isomer)	318
15	Triisopropyl phosphate	224
16	Diisopropyl phosphoric ethyl isopropyl phosphoric anhydride	332



LC-ESI-MS total-ion-current (700 to 70 Da) chromatogram obtained for the munitions-grade tabun sample with a sampling cone voltage of 30 volts. Gradient program: 1% to 75%B in 45 minutes.

## LC-ESI-MS sensitivity

The time-of-flight mass spectrometer instrument offers a significant improvement in sensitivity over the older DRES Autospec-Q instrument for the acquisition of complete electrospray mass spectra. The sensitivity of the LC-ESI-MS method was estimated with triethyl phosphate (Figure 20), a compound that does not readily hydrolyse in water. A complete interpretable mass spectrum (sampling cone voltage: 20 volts) was obtained for 50 pg of TEP during LC-ESI-MS analysis. The S/N ratio for the  $m/z$  183,  $(M+H)^+$ , reconstructed ion current chromatogram was approximately 25:1 (24).

The relative sensitivities of LC-ESI-MS (time-of-flight) and GC-MS (quadrupole) during the acquisition of full mass spectra were investigated with a synthetic sample that failed to produce tabun (26). A 5  $\mu$ L aliquot of the synthetic sample was diluted to 30 mL with water for LC-ESI-MS analysis and to 30 mL of dichloromethane for GC-MS analysis. LC-ESI-MS (5  $\mu$ L injection volume) sample sensitivity for the acquisition of a full mass spectrum was comparable to GC-MS (1  $\mu$ L injection volume) for the compounds detected by both techniques. Similar S/N ratios in the total-ion-current were observed for trace components (low to subnanogram quantities) in the sample using the injection volumes indicated. Twelve sample components were identified by LC-ESI-MS, while only ten were detected by GC-MS. The two sample components that were not detected, ethyl phosphoric tetramethylphosphorodiamidic anhydride (Figure 43) and octamethyltetramidotriphosphoric acid (Figure 59), were not detected by GC-MS due to their polarity. Derivatization and a second analysis would be required for GC-MS analysis of these compounds.

## ESI-MS mass spectral differences

ESI data may be acquired at different sampling cone voltages. At lower sampling cone voltages (e.g., 20 volts) the acquired mass spectra contain molecular mass information in the form of  $MH^+$ ,  $[MH+ACN]^+$  ions and/or protonated dimers, while higher sampling cone voltages (e.g., 30 to 70 volts) will reduce the relative intensity of the  $MH^+$  and enhance the formation of structurally important product ions in the ESI interface. GC-MS data acquired under electron impact ionization generally contain a larger number of structurally significant fragmentation ions, but often these mass spectra do not contain a molecular ion that may be used to establish the molecular mass of the compound. In these cases a second GC-MS analysis under chemical ionization conditions would be required to establish the molecular mass of the identified compound.

Organophosphorous chemical warfare agents and related compounds containing alkoxy substituents generally exhibit ESI-MS product ions due to alkene loss from the alkoxy substituents (and their acetonitrile adducts) at higher sampling cone voltages. Diisopropyl methylphosphonate (Figure 16), an organophosphorus compound identified in a snow sample contaminated with sarin (23), was typical, exhibiting product ions at  $m/z$  139 and  $m/z$  97 due to sequential loss of  $C_3H_6$  from the  $MH^+$  ion.

A common degradation product in VX, bis[2-(diisopropylamino)ethyl] disulfide (Figure 55), was characterized with sampling cone voltages of 20 and 50 volts. Molecular mass was established by the presence of a well resolved  $[\text{MH}_2]^{2+}$  isotopic cluster and a  $\text{MH}^+$  ion at  $m/z$  161 (rounded up to nominal mass) and  $m/z$  321, respectively. Bis[2-(diisopropylamino)ethyl] disulfide, one of a series of compounds of the general form,  $(\text{iPr})_2\text{NC}_2\text{H}_4\text{S}_n\text{C}_2\text{H}_4\text{N}(\text{iPr})_2$  (where  $n=1,2,3$ ), exhibited a number of characteristic product ions due to loss of  $\text{HS}_2\text{C}_2\text{H}_4\text{N}(\text{iPr})_2$ ,  $\text{HSC}_2\text{H}_4\text{N}(\text{iPr})_2$  and  $\text{HN}(\text{iPr})_2$  from the  $\text{MH}^+$  ion at the higher sampling cone voltage.

Hexamethylphosphorotriamide (Figure 13), a product observed in the failed tabun synthetic sample (26) produced only one significant product ion at a higher sampling cone voltage. An ion at  $m/z$  135, due to loss of  $\text{HN}(\text{CH}_3)_2$  from the  $\text{MH}^+$  ion, was observed for this and several other similar amides.

The hydrolysis products of munitions grade mustard included a number of longer chain diols. Diol ESI-MS data contained both molecular ion and product ion content, enabling structural identification of these hydrolysis products. The ESI-MS data obtained for 6-oxa-3,9-dithia-1,11-undecanediol (Figure 37), the hydrolysis product of bis[(2-chloroethylthio)ethyl]ether (T), was typical. The mass spectrum contained a significant  $\text{MH}^+$  ion at  $m/z$  227 and a  $\text{M}_2\text{H}^+$  ion at  $m/z$  453. 6-Oxa-3,9-dithia-1,11-undecanediol exhibited a product ion due to loss of  $\text{H}_2\text{O}$  at  $m/z$  209 along with product ions at  $m/z$  181,  $m/z$  149 and  $m/z$  105 due to  $[\text{MH}-\text{H}_2\text{O}-\text{C}_2\text{H}_4]^+$ ,  $[\text{MH}-\text{H}_2\text{O}-\text{SC}_2\text{H}_4]^+$  and  $[\text{MH}-\text{H}_2\text{O}-\text{SC}_2\text{H}_4-\text{OC}_2\text{H}_4]^+$ , respectively, that could be used to establish S and O positioning.

During the course of investigations a number of common lower mass product ions were observed for the organophosphorus and organosulfur compounds analysed. The presence of these ions, listed in Table 3, additional product ions and the  $\text{MH}^+$  ion often aided in the identification of unknowns. For example, methylphosphonofluoridates such as cyclohexyl methylphosphonofluoridate (Figure 14) tend to exhibit a significant product ion at  $m/z$  99 due to alkene loss from the alkoxy substituent, while phosphates, such as triethyl phosphate (Figure 20) will also exhibit an ion of the same nominal mass (with less relative intensity) after sequential loss of the alkenes associated with the alkoxy substituents. Similarly the presence of a significant ion at  $m/z$  97 often indicates the presence of an alkyl methylphosphonic acid such as pinacolyl methylphosphonic acid (Figure 15) or a dialkyl methylphosphonate such as diethyl methylphosphonate (Figure 7).

Compounds related to VX, such as S-[2-(diisopropylamino) ethyl]methylphosphonothiolate (Figure 39), generally exhibit a common product ion at  $m/z$  128 due to the presence of  $-\text{CH}_2\text{CH}_2\text{N}(\text{CH}(\text{CH}_3)_2)_2$ , while thiodiglycol (Figure 2) and the hydrolysis products of longer chain sulfur sulfur vesicants, such as 3,6,9-trithia-1,11-undecanediol (Figure 41) usually produce a significant protonated thioxane ion at  $m/z$  105. Tabun related compounds, such as diisopropyl dimethylphosphoramidate (Figure 32), containing a single dimethyl amine substituent (and alkoxy substitution) produce lower mass product ions at  $m/z$  126, while an ion at  $m/z$  135 generally indicates an organophosphorus compound with two dimethyl amine substituents.

**Table 3.** Lower mass product ions observed during chemical warfare agent analysis by LC-ESI-MS.

m/z	Possible structures	m/z	Possible structures
97	$\begin{array}{c} \text{OH} \\   \\ \text{CH}_3-\text{P}^+-\text{OH} \\   \\ \text{OH} \end{array}$	127	$\begin{array}{c} \text{OH} \\   \\ \text{CH}_3\text{CH}_2\text{O}-\text{P}^+-\text{OH} \\   \\ \text{OH} \end{array}$
99	$\begin{array}{c} \text{OH} \\   \\ \text{HO}-\text{P}^+-\text{OH} \\   \\ \text{OH} \end{array} \quad \text{or} \quad \begin{array}{c} \text{OH} \\   \\ \text{F}-\text{P}^+-\text{OH} \\   \\ \text{CH}_3 \end{array}$	128	$\begin{array}{c} (\text{CH}_3)_2\text{CH}-\text{N}^+-\text{CH}_2 \\   \quad \quad   \\ (\text{CH}_3)_2\text{CH}-\text{CH}_2 \end{array}$
101	$\begin{array}{c} \text{OH} \\   \\ \text{F}-\text{P}^+-\text{OH} \\   \\ \text{OH} \end{array}$	135	$\begin{array}{c} \text{O} \\    \\ (\text{CH}_3)_2\text{N}-\text{P}^+- \\   \\ \text{N}(\text{CH}_3)_2 \end{array}$
105	$\begin{array}{c} \text{H} \\   \\ \text{S}^+ \\ / \quad \backslash \\ \text{H}_2\text{C} \quad \text{CH}_2 \\   \quad \quad   \\ \text{H}_2\text{C} \quad \text{CH}_2 \\ \backslash \quad / \\ \text{O} \end{array}$	153	$\begin{array}{c} \text{OH} \\   \\ (\text{CH}_3)_2\text{N}-\text{P}^+-\text{OH} \\   \\ \text{N}(\text{CH}_3)_2 \end{array}$
111	$\begin{array}{c} \text{OH} \\   \\ \text{CH}_3\text{CH}_2-\text{P}^+-\text{OH} \\   \\ \text{OH} \end{array}$	154	$\begin{array}{c} \text{OH} \\   \\ \text{CH}_3\text{CH}_2\text{O}-\text{P}^+-\text{OH} \\   \\ \text{N}(\text{CH}_3)_2 \end{array}$
126	$\begin{array}{c} \text{OH} \\   \\ (\text{CH}_3)_2\text{N}-\text{P}^+-\text{OH} \\   \\ \text{OH} \end{array}$	155	$\begin{array}{c} \text{OH} \\   \\ \text{CH}_3\text{CH}_2\text{O}-\text{P}^+-\text{OH} \\   \\ \text{OCH}_2\text{CH}_3 \end{array}$

Accurate mass measurement during LC-ESI-MS analysis may also aid in the assignment of elemental composition of both  $\text{MH}^+$  and product ions. High resolution measurements (5000 resolution with a 50% valley) were made on all the compounds identified in the contaminated snow sample (23). Data acquired under these conditions compared favorably with theoretical data, with errors typically being less than 0.005 Da. Table 4 contains the  $\text{MH}^+$  and selected product ion data for several of the significant sample components identified during LC-ESI-MS of this sample.

**Table 4.** Typical high resolution ESI-MS data acquired for compounds identified in the contaminated snow sample. (Sampling cone voltage: 20 volts) (23).

Compound Name	Ion	Observed Mass (Da) (mean $\pm$ SD; n = 5)	Theoretical Mass (Da)	Average Error (Da)
Methylphosphonic acid	MH <sup>+</sup>	97.0099 $\pm$ 0.0044	97.0055	0.0044
Isopropyl methylphosphonic acid	MH <sup>+</sup>	139.0455 $\pm$ 0.0019	139.0524	0.0069
	[MH-C <sub>3</sub> H <sub>6</sub> ] <sup>+</sup>	97.0071 $\pm$ 0.0015	97.0055	0.0016
Sarin (Isopropyl methylphosphonofluoridate)	MH <sup>+</sup>	141.0444 $\pm$ 0.0021	141.0481	0.0037
	[MH-C <sub>3</sub> H <sub>6</sub> ] <sup>+</sup>	99.0036 $\pm$ 0.0014	99.0011	0.0025
Diisopropyl methylphosphonate	MH <sup>+</sup>	181.1004 $\pm$ 0.0028	181.0994	0.0010
	[MH-C <sub>3</sub> H <sub>6</sub> ] <sup>+</sup>	139.0547 $\pm$ 0.0019	139.0524	0.0023
Tributylamine	MH <sup>+</sup>	186.2221 $\pm$ 0.0026	186.2222	0.0001
Diisopropyl phosphorofluoridate <sup>2</sup>	MH <sup>+</sup>	185.0747 $\pm$ 0.0003	185.0743	0.0004
	[MH-C <sub>3</sub> H <sub>6</sub> ] <sup>+</sup>	143.0285 $\pm$ 0.0010	143.0273	0.0012
Triisopropyl phosphate	MH <sup>+</sup>	225.1246 $\pm$ 0.0027	225.1256	0.0010
	[MH-C <sub>3</sub> H <sub>6</sub> ] <sup>+</sup>	183.0814 $\pm$ 0.0005	183.0786	0.0028


Cyclohexyl methylphosphonofluoridate (GF) (Figure 14), pinacolyl methylphosphonic acid (Figure 15), diisopropyl methylphosphonate (Figure 16) and ethyl tetramethylphosphorodiamidate (Figure 17), four compounds with the same nominal mass (molecular mass of 180), produce a common MH<sup>+</sup> ion at m/z 181, but may be readily differentiated on the basis of their product ions. Cyclohexyl methylphosphonofluoridate exhibited the characteristic methylphosphonofluoridate ion at m/z 99, due to loss of C<sub>6</sub>H<sub>10</sub>, (and its acetonitrile adduct) at m/z 140. Both pinacolyl methylphosphonic acid and diisopropyl methylphosphonate, two compounds with identical elemental composition, exhibited the characteristic m/z 97 product ion (and its acetonitrile adduct) indicating the presence of an alkyl methylphosphonic acid or dialkyl methylphosphonate. The m/z 97 was due to loss of C<sub>6</sub>H<sub>12</sub> from the MH<sup>+</sup> ion of pinacolyl methylphosphonic acid, while diisopropyl methylphosphonate exhibited sequential losses of C<sub>3</sub>H<sub>6</sub>, resulting in an additional product ion at m/z 139. The final compound with a molecular mass of 180 Da, ethyl tetramethylphosphorodiamidate, exhibited several characteristic product ions, the most significant being at m/z 153, due to loss of C<sub>2</sub>H<sub>4</sub> from the MH<sup>+</sup> ion.

Two additional compounds with identical elemental composition, bis(ethyl dimethylphosphoramidic) anhydride (Figure 50) and diethyl phosphoric tetramethylphosphorodiamidic anhydride (Figure 51), exhibited identical mass spectra dominated by a MH<sup>+</sup> ion at m/z 289 under low sampling cone voltage conditions. Both anhydrides contained two ethoxy and two dimethylamine substituents based on the product ions observed. Differentiation of the two possible isomers was possible on the basis of the characteristic lower mass ions observed with a sampling cone voltage of 40 volts. Bis(ethyl dimethylphosphoramidic) anhydride exhibited ions at m/z 126 and m/z 154 indicating the

presence of both ethoxy and dimethyl amine substitution at each phosphorus atom, while the product ions at  $m/z$  135 and  $m/z$  153 for diethyl phosphoric tetramethylphosphorodiamidic anhydride indicated that one phosphorus contained both dimethylamine substituents. Data were obtained under high resolution conditions and Table 5 lists the elemental assignments for the ions observed during accurate mass measurement of both anhydrides (26).


**Table 5.** High resolution ESI-MS data for two compounds identified in a failed tabun synthesis sample (Sampling cone voltage: 40 volts) (26).

Compound Name	Ion	Observed Mass (Da) (mean $\pm$ SD) <sup>2</sup>	Theoretical Mass (Da)	Average Error (Da)
Diethyl phosphoric tetramethylphosphorodiamidic anhydride	MH <sup>+</sup>	289.1049 $\pm$ 0.0011	289.1082	0.0033
	[MH-C <sub>2</sub> H <sub>4</sub> ] <sup>+</sup>	261.0790 $\pm$ 0.0006	261.0769	0.0021
	[MH-HNMe <sub>2</sub> ] <sup>+</sup>	244.0554 $\pm$ 0.0020	244.0504	0.0050
	[MH-(C <sub>2</sub> H <sub>4</sub> ) <sub>2</sub> ] <sup>+</sup>	233.0489 $\pm$ 0.0007	233.0456	0.0033
	[(Me <sub>2</sub> N) <sub>2</sub> P(OH) <sub>2</sub> ] <sup>+</sup>	153.0826 $\pm$ 0.0012	153.0793	0.0033
	[(Me <sub>2</sub> N) <sub>2</sub> PO] <sup>+</sup>	135.0705 $\pm$ 0.0009	135.0687	0.0018
Bis(ethyl dimethylphosphoramidic) anhydride	MH <sup>+</sup>	289.1064 $\pm$ 0.0011	289.1082	0.0018
	[MH-C <sub>2</sub> H <sub>4</sub> ] <sup>+</sup>	261.0753 $\pm$ 0.0015	261.0769	0.0016
	[MH-HNMe <sub>2</sub> ] <sup>+</sup>	244.0517 $\pm$ 0.0021	244.0504	0.0013
	[MH-HNMe <sub>2</sub> -C <sub>2</sub> H <sub>4</sub> ] <sup>+</sup>	229.0112 $\pm$ 0.0024	229.0143	0.0031
	[(Me <sub>2</sub> N)(EtO)P(OH) <sub>2</sub> ] <sup>+</sup>	154.0648 $\pm$ 0.0018	154.0633	0.0015
	[(Me <sub>2</sub> N)P(OH) <sub>3</sub> ] <sup>+</sup>	126.0374 $\pm$ 0.0019	126.0320	0.0054



**Figure 1. Diisopropylamine (101)**


CAS: 108-18-9  
 Formula:  $C_6H_{15}N$   
 Monoisotopic MW: 101.1204  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™



**Figure 2. 2,2'-Thiodiethanol (Thiodiglycol) (122)**



CAS: 111-48-8  
 Formula:  $C_4H_{10}O_2S$   
 Monoisotopic MW: 122.0402  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™






**Figure 3. Ethyl methylphosphonic acid (124)**

CAS: 1832-53-7  
 Formula:  $\text{C}_3\text{H}_9\text{O}_3\text{P}$   
 Monoisotopic MW: 124.0289  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™




**Figure 4. Isopropyl methylphosphonic acid (138)**

CAS: 1832-54-8  
 Formula: C<sub>4</sub>H<sub>11</sub>O<sub>3</sub>P  
 Monoisotopic MW: 138.0446  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™




**Figure 5. Ethyl ethylphosphonic acid (138)**

CAS:  
 Formula:  $C_4H_{11}O_3P$   
 Monoisotopic MW: 138.0446  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™




**Figure 6. Isopropyl methylphosphonofluoridate (Sarin, GB) (140)**

CAS: 107-44-8  
 Formula:  $\text{C}_4\text{H}_{10}\text{FO}_2\text{P}$   
 Monoisotopic MW: 140.0403  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™





**Figure 7. Diethyl methylphosphonate (152)**

CAS: 683-08-9  
 Formula:  $C_5H_{13}O_3P$   
 Monoisotopic MW: 152.0602  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 8. Isopropyl methyl methylphosphonate (152)**

CAS: 690-64-2  
 Formula:  $C_5H_{13}O_3P$   
 Monoisotopic MW: 152.0602  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™



**Figure 9. Ethyl dimethylphosphoramidic acid (153)**



CAS: 2632-86-2  
 Formula:  $C_4H_{12}NO_3P$   
 Monoisotopic MW: 153.0555  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™



**Figure 10. Tetramethyl phosphorodiamidic cyanide (161)**


CAS: 14445-60-4  
 Formula:  $C_5H_{12}N_3OP$   
 Monoisotopic MW: 161.0718  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™







**Figure 11. O-Ethyl N,N-dimethylphosphoroamidocyanidate (Tabun, GA) (162)**

CAS: 77-81-6  
 Formula: C<sub>5</sub>H<sub>11</sub>N<sub>2</sub>O<sub>2</sub>P  
 Monoisotopic MW: 162.00558  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™





**Figure 12. Cyclohexyl methylphosphonic acid (178)**

CAS: 1932-60-1  
 Formula:  $C_7H_{15}O_3P$   
 Monoisotopic MW: 178.0759  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 13. Hexamethylphosphorotriamide (179)**

CAS: 680-31-9  
 Formula:  $C_6H_{18}N_3OP$   
 Monoisotopic MW: 179.1188  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 14. Cyclohexyl methylphosphonofluoridate (GF) (180)**

CAS: 329-99-7  
 Formula:  $C_7H_{14}FO_2P$   
 Monoisotopic MW: 180.0715  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 15. Pinacolyl methylphosphonic acid (180)**

CAS: 616-52-4  
 Formula: C<sub>7</sub>H<sub>17</sub>O<sub>3</sub>P  
 Monoisotopic MW: 180.0915  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™





**Figure 16. Diisopropyl methylphosphonate (180)**

CAS: 1445-75-6  
 Formula:  $C_7H_{17}O_3P$   
 Monoisotopic MW: 180.0915  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™



**Figure 17. Ethyl tetramethylphosphorodiamidate (180)**



CAS: 2404-65-1  
 Formula:  $C_6H_{17}N_2O_2P$   
 Monoisotopic MW: 180.1028  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™



**Figure 18. Diethyl dimethylphosphoramidate (181)**



CAS: 2404-03-7  
 Formula:  $C_6H_{16}NO_3P$   
 Monoisotopic MW: 181.0868  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™






**Figure 19. 3,6-Dithia-1,8-octanediol (182)**

CAS: 5244-34-8  
 Formula: C<sub>6</sub>H<sub>14</sub>O<sub>2</sub>S<sub>2</sub>  
 Monoisotopic MW: 182.0435  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™





**Figure 20. Triethyl phosphate (182)**

CAS: 78-40-0  
 Formula:  $C_6H_{15}O_4P$   
 Monoisotopic MW: 182.0708  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 21. Pinacolyl methylphosphonofluoridate (Soman, GD) (182)**

CAS: 96-64-0  
 Formula:  $\text{C}_7\text{H}_{16}\text{FO}_2\text{P}$   
 Monoisotopic MW: 182.0872  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™




**Figure 22. 8-Chloro-6-oxa-3-thia-1-octanol (184)**

CAS: Unknown  
 Formula:  $C_6H_{13}ClO_2S$   
 Monoisotopic MW: 184.0325  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™




**Figure 23. Diisopropyl phosphorofluoridate (DFP) (184)**

CAS: 55-91-4  
 Formula:  $\text{C}_6\text{H}_{14}\text{FO}_3\text{P}$   
 Monoisotopic MW: 184.0665  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™




**Figure 24. Tributylamine (185)**

CAS: 102-82-9  
 Formula:  $\text{C}_{12}\text{H}_{27}\text{N}$   
 Monoisotopic MW: 185.2144  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™



**Figure 25. 2-(Diisopropylamino)ethyl ethyl sulfide (189)**



CAS:  
 Formula:  $C_{10}H_{23}NS$   
 Monoisotopic MW: 189.1551  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™



**Figure 26. Diisopropyl ethylphosphonate (194)**



CAS: 1067-69-2  
 Formula:  $C_8H_{19}O_3P$   
 Monoisotopic MW: 194.1072  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™







**Figure 27. Ethyl isopropyl dimethylphosphoramidate (195)**

CAS: 99520-56-6  
 Formula:  $\text{C}_7\text{H}_{18}\text{NO}_3\text{P}$   
 Monoisotopic MW: 195.1024  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™





**Figure 28. Diethyl isopropyl phosphate (196)**

CAS: 2736-99-4  
 Formula: C<sub>7</sub>H<sub>17</sub>O<sub>4</sub>P  
 Monoisotopic MW: 196.0864  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™




**Figure 29. Pentyl dibutylamine (199)**

CAS: 6703-22-6  
 Formula:  $C_{13}H_{29}N$   
 Monoisotopic MW: 199.2300  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™




**Figure 30. 2,2',2''-Trichlorotriethylamine (HN3) (203)**

CAS: 555-77-1  
 Formula:  $C_6H_{12}Cl_3N$   
 Monoisotopic MW: 203.0035  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 31. 1,7-Dioxo-4,10-dithiacyclodecane (208)**

CAS: 294-95-1  
Formula:  $C_8H_{16}O_2S_2$   
Monoisotopic MW: 208.0592  
Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 32. Diisopropyl dimethylphosphoramidate (209)**

CAS: 2404-04-8  
 Formula: C<sub>8</sub>H<sub>20</sub>NO<sub>3</sub>P  
 Monoisotopic MW: 209.1181  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™



**Figure 33. Diisopropyl ethyl phosphate (210)**



CAS: 2695-64-6  
 Formula:  $\text{C}_8\text{H}_{19}\text{O}_4\text{P}$   
 Monoisotopic MW: 210.1021  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™



**Figure 34. Diisopropyl chloromethylphosphonate (214)**



CAS: 6954-83-2  
 Formula: C<sub>7</sub>H<sub>16</sub>ClO<sub>3</sub>P  
 Monoisotopic MW: 214.0526  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™







**Figure 35. Triisopropyl phosphate (224)**

CAS: 513-08-6  
 Formula:  $\text{C}_9\text{H}_{21}\text{O}_4\text{P}$   
 Monoisotopic MW: 224.1177  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™





**Figure 36. *N,N'*-Dicyclohexylurea (224)**

CAS: 2387-23-7  
 Formula:  $C_{13}H_{24}N_2O$   
 Monoisotopic MW: 224.1889  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 37. 6-Oxa-3,9-dithia-1,11-undecanediol (226)**

CAS: 7426-02-0  
 Formula:  $C_8H_{18}O_3S_2$   
 Monoisotopic MW: 226.0697  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™




**Figure 38. Diethyl dimethylpyrophosphonate (230)**


CAS: 32288-17-8  
 Formula:  $\text{C}_6\text{H}_{16}\text{O}_5\text{P}_2$   
 Monoisotopic MW: 230.0473  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™



**Figure 39. *S*-[2-(Diisopropylamino) ethyl]methylphosphonothiolate (239)**



CAS: 73207-98-4  
 Formula:  $C_9H_{22}NO_2PS$   
 Monoisotopic MW: 239.1109  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 41. 3,6,9-Trithia-1,11-undecanediol (242)**



CAS: 14440-77-8  
 Formula:  $C_8H_{18}O_2S_3$   
 Monoisotopic MW: 242.0469  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™



**Figure 42. Diisopropyl dichloromethylphosphonate (248)**



CAS: 55696-11-2  
 Formula:  $\text{C}_7\text{H}_{15}\text{Cl}_2\text{O}_3\text{P}$   
 Monoisotopic MW: 248.0136  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™







**Figure 43. Ethyl phosphoric tetramethylphosphorodiamidic anhydride (260)**

CAS:  
Formula: C6H18N2O5P2  
Monoisotopic MW: 260.0691  
Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™




**Figure 44.** Diethyl phosphoric ethyl hydrogen phosphoric anhydride (262)

CAS: 20680-55-1  
 Formula: C<sub>6</sub>H<sub>16</sub>O<sub>7</sub>P<sub>2</sub>  
 Monoisotopic MW: 262.0371  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™





**Figure 45. *O*-Ethyl *S*-[2-(diisopropylamino)ethyl] methylphosphonothiolate (VX) (267)**

CAS: 50782-69-9  
 Formula:  $C_{11}H_{26}NO_2PS$   
 Monoisotopic MW: 267.1422  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 46. 6-Oxa-3,9,12-trithia-1,4-tetradecanediol (286)**




**Figure 47. Bis(tetramethylphosphorodiamidic) anhydride (286)**

CAS: 152-16-9  
 Formula:  $\text{C}_8\text{H}_{24}\text{N}_4\text{O}_3\text{P}_2$   
 Monoisotopic MW: 286.1324  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™





**Figure 48. Ethyl dimethylphosphoramidic tetramethylphosphoamidic anhydride (287)**

CAS: 29506-48-7  
 Formula: C<sub>8</sub>H<sub>23</sub>N<sub>3</sub>O<sub>4</sub>P<sub>2</sub>  
 Monoisotopic MW: 287.1164  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™



**Figure 49. Bis[2-(diisopropylamino)ethyl] sulfide (288)**



CAS: 110501-56-9  
 Formula:  $\text{C}_{16}\text{H}_{36}\text{N}_2\text{S}$   
 Monoisotopic MW: 288.2599  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™



**Figure 50. Bis(ethyl dimethylphosphoramidic) anhydride (288)**



CAS: 28616-48-0  
 Formula: C<sub>8</sub>H<sub>22</sub>N<sub>2</sub>O<sub>5</sub>P<sub>2</sub>  
 Monoisotopic MW: 288.1004  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™






**Figure 51. Diethyl phosphoric tetramethylphosphorodiamidic anhydride(288)**

CAS: 1474-80-2  
 Formula:  $\text{C}_8\text{H}_{22}\text{N}_2\text{O}_5\text{P}_2$   
 Monoisotopic MW: 288.1004  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™




**Figure 52. Diethyl phosphoric ethyl dimethylphosphoramidic anhydride (289)**

CAS: 92403-54-8  
 Formula: C<sub>8</sub>H<sub>21</sub>NO<sub>6</sub>P<sub>2</sub>  
 Monoisotopic MW: 289.0844  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™





**Figure 53. Bis(diethyl phosphoric) anhydride (TEPP) (290)**

CAS: 107-49-3  
 Formula:  $\text{C}_8\text{H}_{20}\text{O}_7\text{P}_2$   
 Monoisotopic MW: 290.0684  
 Instrumental: LC-ESI-MS (ACN/H<sub>2</sub>O) using a Micromass LCT™





**Figure 54. Diethyl phosphoric ethyl isopropyl phosphoric anhydride (304)**

CAS: 120551-44-2  
 Formula:  $\text{C}_9\text{H}_{22}\text{O}_7\text{P}_2$   
 Monoisotopic MW: 304.0841  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™




**Figure 55. Bis[2-diisopropylamino)ethyl] disulfide (320)**

CAS: 65332-44-7  
 Formula:  $C_{16}H_{36}N_2S_2$   
 Monoisotopic MW: 320.2320  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™





**Figure 56. O-Ethyl-S-[5-(diisopropylamino)-3-thiapentyl] methylphosphonothiolate (327)**

CAS: 110501-57-0  
 Formula:  $\text{C}_{13}\text{H}_{30}\text{NO}_2\text{PS}_2$   
 Monoisotopic MW: 327.1456  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™



**Figure 57. Diisopropyl phosphoric ethyl isopropyl phosphoric anhydride (332)**


CAS: 63704-58-5  
 Formula:  $\text{C}_{11}\text{H}_{26}\text{O}_7\text{P}_2$   
 Monoisotopic MW: 332.1154  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™



**Figure 58. Bis(diisopropylphosphoric) anhydride (346)**


CAS: 5836-28-2  
 Formula:  $\text{C}_{12}\text{H}_{28}\text{O}_7\text{P}_2$   
 Monoisotopic MW: 346.1310  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™





**Figure 59. Octamethyltetramidotriphosphoric acid (or isomer) (366)**

CAS:  
 Formula:  $\text{C}_8\text{H}_{25}\text{N}_4\text{O}_6\text{P}_3$   
 Monoisotopic MW: 366.0987  
 Instrumental: LC-ESI-MS (ACN/ $\text{H}_2\text{O}$ ) using a Micromass LCT™



**Figure 60. O-Ethyl-S[8-(diisopropylamino)-3,6-dithiaoctyl]methylphosphonothiolate (387)**

CAS: 110501-60-5  
 Formula:  $C_{15}H_{34}NO_2PS_3$   
 Monoisotopic MW: 387.1489  
 Instrumental: LC-ESI-MS (ACN/ $H_2O$ ) using a Micromass LCT™

## Conclusions

---

LC-ESI-MS provides an advantage over GC-MS for the analysis of aqueous samples and extracts, as both the chemical warfare agents and their hydrolysis products may be analysed at the nanogram to picogram level during a single analysis without the need for additional sample handling or derivatization. Over the past several years DRES developed LC-ESI-MS methods for the analysis of a variety of contaminated samples including; hydrolysed munitions grade mustard, degraded O-ethyl S-[2-(diisopropylamino)ethyl] methylphosphonothiolate (VX), snow contaminated with sarin (GB), soil contaminated with sarin (GB) and soman (GD), munitions grade tabun (GA), failed synthetic products and synthetic standards. During the course of these investigations DRES acquired ESI-MS data for sixty chemical warfare agents, degradation products and related compounds. Data were obtained with a resolution of 5000 (50% valley definition) in the continuum mode at several sampling cone voltages. At lower sampling cone voltages the mass spectra were dominated by protonated, sodiated and protonated acetonitrile adducts and/or their dimers that could be used to confirm the molecular mass of each compound. Structural information was obtained by inducing product ion formation in the ESI interface at higher sampling cone voltages.

Representative data obtained at both lower and higher sampling cone voltages were selected for entry into the DRES ESI-MS Database as part of a collaboration between Canada (DRES) and The Netherlands (TNO Prins Maurits Laboratory). The reported data should prove valuable for the identification of chemical warfare agents, degradation products and related compounds samples collected by the military or in support of Chemical Weapons Convention inspections.

## References

---

1. Kientz, Ch. E. (1998). Chromatography and mass spectrometry of chemical warfare agents, toxins and related compounds: state of the art and future prospects. *J. Chromatogr. A*, 814, 1-23.
2. Witkiewicz, Z., Mazurek, M. and Szulc, J. (1990). Chromatographic analysis of chemical warfare agents. *J. Chromatogr.*, 503, 293-357.
3. Slobodnik, J., van Baar, B. L. M. and Brinkman, U. A. Th. (1995). Column liquid chromatography-mass spectrometry: Selected techniques in environmental applications for polar pesticides and related compounds. *J. Chromatogr. A*, 703, 81-121.
4. Bell, A. J., Despeyroux, D., Murrell, J. and Watts, P. (1997). Fragmentation and reactions of organophosphate ions produced by electrospray ionization. *Int. J. Mass Spectrom. Ion. Proc.*, 165/166, 533-550.
5. Kostiainen, R., Bruins, A. P. and Hakkinen, V. M. A. (1993). Identification of degradation products of some chemical warfare agents by capillary electrophoresis-ion spray mass spectrometry. *J. Chromatogr.*, 634, 113-118.
6. D'Agostino, P. A., Provost, L. R. and Hancock, J. R. (1994). Mass spectrometric identification of chemical warfare agent precursor and degradation products scheduled by the chemical weapons convention, In *Proceedings of the 42 nd Annual Conference on Mass Spectrometry and Allied Topics*, Chicago, Illinois, 275-276.
7. Borrett, V. T., Colton, R. and Traeger, J. C. (1995). The electrospray mass spectra of phosphonic acid, methyl phosphonic acid and its alkyl esters, and their complexes with alkali and alkali earth metal ions. *Eur. Mass Spectrom.*, 1, 131-140.
8. Borrett, V. T., Mathews, R. J., Colton, R. and Traeger, J. C. (1996). Verification of the United Nations chemical weapons convention: The application of electrospray mass spectrometry. *Rapid Commun. Mass Spectrom.*, 10, 114-118.
9. Black, R. M. and Read, R. W. (1997). Application of liquid chromatography-atmospheric pressure chemical ionization mass spectrometry, and tandem mass spectrometry, to the analysis and identification of degradation products of chemical warfare agents. *J. Chromatogr. A*, 759, 79-92.
10. Black, R. M. and Read, R. W. (1998). Analysis of degradation products of organophosphorus chemical warfare agents and related compounds by liquid chromatography-mass spectrometry using electrospray and atmospheric pressure chemical ionization. *J. Chromatogr. A*, 794, 233-244.

11. D'Agostino, P. A., Provost, L. R. and Hancock, J. R. (1998). Packed Capillary Column Electrospray Mass Spectrometry and Tandem Mass Spectrometry of Hydrolysed HT and HQ, Suffield Report No 691, 23 pages.
12. D'Agostino, P. A., Provost, L. R. and Hancock, J. R. (1998). Analysis of mustard hydrolysis products by packed capillary liquid chromatography-electrospray mass spectrometry. *J. Chromatogr. A*, 808, 177-184.
13. Mercier, J.-P., Chaimbault, P., Morin, Ph., Dreux, M. and Tambute, A. (1998). Identification of phosphonic acids by capillary electrophoresis-ion spray mass spectrometry. *J. Chromatogr. A*, 825, 71-80.
14. D'Agostino, P. A., Hancock, J. R. and Provost, L. R. (1999). Packed Capillary LC-ESI-MS Analysis of O-Ethyl S-[2-(diisopropylamino)ethyl] Methylphosphonothiolate (VX), Suffield Report No 706, 28 pages.
15. D'Agostino, P. A., Hancock, J. R. and Provost, L. R. (1999). Analysis of O-ethyl S-[2-(diisopropylamino)ethyl] methylphosphonothiolate (VX) and its degradation products by packed capillary liquid chromatography-electrospray mass spectrometry. *J. Chromatogr. A*, 837, 93-105.
16. D'Agostino, P. A., Hancock, J. R. and Provost, L. R. (1999). Analysis of GA, GB, GD and GF in Aqueous Samples by Packed Capillary LC-ESI-MS, DRES TM 1999-047, 16 pages.
17. D'Agostino, P. A., Hancock, J. R. and Provost, L. R. (1999). Packed capillary liquid chromatography-electrospray mass spectrometry analysis of organophosphorus chemical warfare agents. *J. Chromatogr. A*, 840, 289-294.
18. Read, R. W. and Black, R. M. (1999). Rapid screening procedures for the hydrolysis products of chemical warfare agents using positive and negative ion liquid chromatography-mass spectrometry and atmospheric pressure chemical ionization. *J. Chromatogr. A*, 862, 169-177.
19. Hooijschuur, E. W., Kientz, C. E., Hulst, A. G. and Brinkman, U. A. Th. (2000). Determination of hydrolysis products of sulfur mustard by reversed-phase microcolumn liquid chromatography coupled on-line with sulfur flame photometric detection and electrospray ionization mass spectrometry using large-volume injections and peak compression. *Anal. Chem.* 72, 1199-1206.
20. D'Agostino, P. A., Hancock, J. R. and Provost, L. R. (2000). Analysis of Tabun and Related Compounds by Packed Capillary Liquid Chromatography-Electrospray Mass Spectrometry (LC-ESI-MS), DRES TM 2000-004, 25 pages.
21. D'Agostino, P. A., Hancock, J. R. and Provost, L. R. (2000). Determination of Organophosphorus Chemical Warfare Agents and their Hydrolysis Products in Soil by Packed Capillary Liquid Chromatography-Electrospray Mass Spectrometry, DRES TR 2000-075, 19 pages.

22. D'Agostino, P. A., Hancock, J. R. and Provost, L. R. (2001). Determination of sarin, soman and their hydrolysis products in soil by packed capillary liquid chromatography-electrospray mass spectrometry. *J. Chromatogr. A*, 912, 291-299.
23. D'Agostino, P. A., Chenier, C. L. and Hancock, J. R. (2001). Identification of sarin and related compounds in snow by packed capillary liquid chromatography-electrospray mass spectrometry. DRES TM 2001-044, 24 pages, UNCLASSIFIED.
24. D'Agostino, P. A., Hancock, J. R. and Provost, L. R. (2001). Electrospray mass spectrometry of chemical warfare agents. In Gelpi E. (Ed), *Advances in Mass Spectrometry*, 15, 297-316. John Wiley and Sons Ltd., New York.
25. van Baar, B. L. M., Hulst, A. G. and Wils, E. R. J. (2001). Electrospray mass spectrometry and tandem mass spectrometry of some alkyl methylphosphonates and ethylphosphonates. TNO Report PML 2000-A71, 20 pages.
26. D'Agostino, P. A. (unpublished results).

UNCLASSIFIED

DOCUMENT CONTROL DATA		
(Security classification of title, body of abstract and indexing annotation must be entered when the overall document is classified)		
1. ORIGINATOR (the name and address of the organization preparing the document. Organizations for who the document was prepared, e.g. Establishment sponsoring a contractor's report, or tasking agency, are entered in Section 8.)  DRE Suffield	2. SECURITY CLASSIFICATION (overall security classification of the document, including special warning terms if applicable)  UNCLASSIFIED	
3. TITLE (the complete document title as indicated on the title page. Its classification should be indicated by the appropriate abbreviation (S, C or U) in parentheses after the title).  Electrospray Mass Spectrometry of Chemical Warfare Agents, Degradation Products and Related Compounds (U)		
4. AUTHORS (Last name, first name, middle initial. If military, show rank, e.g. Doe, Maj. John E.)  D'Agostino, Paul A., Chenier, Claude L. and Hancock, James R.		
5. DATE OF PUBLICATION (month and year of publication of document)  March 2002	6a. NO. OF PAGES (total containing information, include Annexes, Appendices, etc) 84	6b. NO. OF REFS (total cited in document) 26
7. DESCRIPTIVE NOTES (the category of the document, e.g. technical report, technical note or memorandum. If appropriate, enter the type of report, e.g. interim, progress, summary, annual or final. Give the inclusive dates when a specific reporting period is covered.)  Technical Report		
8. SPONSORING ACTIVITY (the name of the department project office or laboratory sponsoring the research and development. Include the address.)		
9a. PROJECT OR GRANT NO. (If appropriate, the applicable research and development project or grant number under which the document was written. Please specify whether project or grant.)	9b. CONTRACT NO. (If appropriate, the applicable number under which the document was written.)	
10a. ORIGINATOR'S DOCUMENT NUMBER (the official document number by which the document is identified by the originating activity. This number must be unique to this document.)  DRES TR 2002-028	10b. OTHER DOCUMENT NOS. (Any other numbers which may be assigned this document either by the originator or by the sponsor.)	
11. DOCUMENT AVAILABILITY (any limitations on further dissemination of the document, other than those imposed by security classification)  <input checked="" type="checkbox"/> Unlimited distribution <input type="checkbox"/> Distribution limited to defence departments and defence contractors; further distribution only as approved <input type="checkbox"/> Distribution limited to defence departments and Canadian defence contractors; further distribution only as approved <input type="checkbox"/> Distribution limited to government departments and agencies; further distribution only as approved <input type="checkbox"/> Distribution limited to defence departments; further distribution only as approved <input type="checkbox"/> Other (please specify):		
12. DOCUMENT ANNOUNCEMENT (any limitation to the bibliographic announcement of this document. This will normally corresponded to the Document Availability (11). However, where further distribution (beyond the audience specified in 11) is possible, a wider announcement audience may be selected).		

UNCLASSIFIED

UNCLASSIFIED

13. **ABSTRACT** (a brief and factual summary of the document. It may also appear elsewhere in the body of the document itself. It is highly desirable that the abstract of classified documents be unclassified. Each paragraph of the abstract shall begin with an indication of the security classification of the information in the paragraph (unless the document itself is unclassified) represented as (S), (C) or (U). It is not necessary to include here abstracts in both official languages unless the text is bilingual).

Over the past several years Defence Research Establishment Suffield (DRES) developed LC-ESI-MS methods for the analysis of a variety of samples of chemical defence interest. These samples included hydrolysed munitions grade mustard, degraded O-ethyl S-[2-(diisopropylamino)ethyl] methylphosphonothiolate (VX), snow contaminated with sarin (GB), soil contaminated with sarin (GB) and soman (GD), munitions grade tabun (GA), failed synthetic products and synthetic standards. During these investigations DRES acquired ESI-MS data for sixty chemical warfare agents, degradation products and related compounds. Data were obtained with a resolution of 5000 (50% valley definition) in the continuum mode at several sampling cone voltages. At lower sampling cone voltages the mass spectra were dominated by protonated, sodiated and protonated acetonitrile adducts and/or their dimers that could be used to confirm the molecular mass of each compound. Structural information was obtained by inducing product ion formation in the ESI interface at higher sampling cone voltages. Representative data obtained at both a lower and higher sampling cone voltage were selected for entry into the DRES ESI-MS Database as part of a collaboration between Canada (DRES) and The Netherlands (TNO Prins Maurits Laboratory).

14. **KEYWORDS, DESCRIPTORS or IDENTIFIERS** (technically meaningful terms or short phrases that characterize a document and could be helpful in cataloguing the document. They should be selected so that no security classification is required. Identifiers, such as equipment model designation, trade name, military project code name, geographic location may also be included. If possible keywords should be selected from a published thesaurus, e.g. Thesaurus of Engineering and Scientific Terms (TEST) and that thesaurus-identified. If it is not possible to select indexing terms which are Unclassified, the classification of each should be indicated as with the title.)

Liquid chromatography

Mass spectrometry

Electrospray

Chemical warfare agent

Degradation product



**Defence R&D Canada**

is the national authority for providing  
Science and Technology (S&T) leadership  
in the advancement and maintenance  
of Canada's defence capabilities.

**R et D pour la défense Canada**

est responsable, au niveau national, pour  
les sciences et la technologie (S et T)  
au service de l'avancement et du maintien des  
capacités de défense du Canada.



[www.drdc-rddc.dnd.ca](http://www.drdc-rddc.dnd.ca)